

aycan mammography workstation powered by aycan workstation





Key Benefits:

- » Supports Breast Tomosynthesis Objects
- » Fast loading times
- » Intuitive, easy to use
- » Save layouts and settings
- » Integrated prior study view
- » Image sharpness with EIZO monitors
- » Live customer service and support

FDA Cleared. Fast and Consistent Reading of Studies.

aycan's mammography workstation with Advanced Hanging Protocols is optimized for mammography workflow, including the evaluation of dynamic breast MRI studies. It offers excellent image loading times and an intuitive user interface that is easy to use. With the Hanging Protocols you can save display layouts and settings in the 2D and tomosynthesis viewer. You can have it automatically display prior studies, which further increases the speed of the daily workflow. The mammography workstation also offers automatic play-back of tomosynthesis loops.

The EIZO GX550 5MP monochrome monitor offers high brightness of 2,000 cd/m2, which brings out subtle differences in similar shades of gray and sharper rendering of details when displaying monochrome digital mammography images. The optional EIZO RX850 8MP color monitor is ideal for viewing all medical images at once including digital mammography, MRI, US, CT, X-ray, etc. The monitors have FDA 510(k) clearance for both mammography and breast tomosynthesis.

aycan mammography workstation includes:

- » aycan mammography workstation software license
- » EIZO Radiforce GX550 or RX850 monitors
- » RadiCS Monitor calibration kit
- » Apple iMac 27" with 32GB Ram
- » Mammography keypad



Contact us: 585-473-1350 . info@aycan.com . www.aycan.com



About aycan workstation

Dedicated to DICOM images, aycan workstation is one of the fastest, most cost effective, FDA cleared, 64-bit image post-processing software for navigation and visualization of multi-slice/multi-dimensional (3D-reconstruction) studies.